



### **Product description**

Order number	80.1301
Product description	Swabs, in the tube, no medium, (LxØ): 175 x 12 mm, stem material: PS, stem length: 131 mm, swab material: viscose, sterile, 100 piece(s)/case

### **Product characteristics**

Closure type	push cap
Design	in the tube
Type of medium	no medium

### **Size**

Diameter	12 mm
Length of product	175 mm
Length including cap	175 mm
Length excluding cap	155 mm
Length including cap and plunger	166 mm
Length of stem	131 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	white
Swab material	viscose
Stem material	Polystyrene (PS)

### **Purity & certification**

Product category	Medical device
CE certified	CE supplier
Purity standard	sterile
Sterilisation	Gamma irradiation
Safety data sheet	no
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	32000
Depth of case	590 mm
Width of case	220 mm
Height of case	300 mm
Case volume	0.0389 cbm
Weight of product	0.0047 kg
Weight of case	5.94 kg
EAN of inner box	4038917086103
EAN case	4038917234184

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

