



Product description

| | |
|---------------------|--|
| Order number | 80.1363 |
| Product description | Swabs, with tube included, Amies gel, (LxØ): 175 x 12 mm, stem material: aluminium, stem length: 133 mm, swab material: viscose, sterile, 1 piece(s)/blister |

Product characteristics

| | |
|----------------|--------------------|
| Closure type | push cap |
| Design | with tube included |
| Type of medium | Amies gel |

Size

| | |
|----------------------------------|--------|
| Diameter | 12 mm |
| Length of product | 175 mm |
| Length including cap | 175 mm |
| Length excluding cap | 155 mm |
| Length including cap and plunger | 164 mm |
| Length of stem | 133 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | orange |
| Swab material | viscose |
| Stem material | aluminium |

Purity & certification

| | |
|-------------------|-------------------|
| Product category | Medical device |
| CE certified | CE supplier |
| Purity standard | sterile |
| Sterilisation | Gamma irradiation |
| Safety data sheet | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | blister |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 11500 |
| Depth of case | 520 mm |
| Width of case | 450 mm |
| Height of case | 220 mm |
| Case volume | 0.0515 cbm |
| Weight of product | 0.0111 kg |
| Weight of case | 6.5 kg |
| EAN of inner box | 4038917086349 |
| EAN case | 4038917004909 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

