



### **Product description**

Order number	92.181.610
Product description	Rotor, for 36 collection tubes, for micro tube Ø 10.8 mm

### **Product characteristics**

Suitable for	micro tube Ø 10.8 mm
Storage capacity	36

### **Purity & certification**

Product category	no medical device   IVD
Batched	yes

### **Packaging**

Minimum order qty.	1
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	40
Case volume	0.0056 cbm
Weight of product	1.3595 kg
Weight of case	1.36 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Rotor, for 36 collection tubes, for micro tube Ø**  
**10.8 mm**

**Page 2**

EAN of inner box                      4038917156233  
EAN case                                 4038917156226

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**