



Product description

Order number 92.181.611
Product description Rotor, for 22 collection tubes

Product characteristics

Storage capacity 22

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1
Piece(s) / inner box 1
Piece (s) / outer case 1
Piece(s) / pallet 40
Case volume 0.0052 cbm
Weight of product 1.3595 kg
Weight of case 1.36 kg
EAN of inner box 4038917156257
EAN case 4038917156240

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

