



### **Product description**

Order number	95.988.004
Product description	Base station, material: PC, yellow, (LxWxH): 128 x 85.5 x 30.4 mm, suitable for work tray 95.987.002 & 0.2 ml PCR tubes, 5 piece(s)/case

### **Product characteristics**

Suitable for	work tray 95.987.002 & 0.2 ml PCR tubes
--------------	---

### **Size**

Width of product	85.5 mm
Height of product	30.4 mm
Length of product	128 mm

### **Material & colours**

Product material	Polycarbonate (PC)
Colour of product	yellow

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                    no

### **Packaging**

Minimum order qty.                    5  
Type of smallest  
subpackaging                            case  
Piece(s) / inner box                    5  
Piece (s) / outer case                   75  
Piece(s) / pallet                        3600  
Depth of case                            598 mm  
Width of case                            398 mm  
Height of case                            145 mm  
Case volume                              0.0345 cbm  
Weight of product                       0.0302 kg  
Weight of case                            2.83 kg  
EAN of inner box                        4038917124195  
EAN case                                  4038917019156

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

