



Product description

| | |
|---------------------|--|
| Order number | 72.735.002 |
| Product description | PCR single tube, 0.5 ml, PCR Performance Tested, transparent, material: PP, with anti-contamination shield, flat cap, suitable for qPCR, 500 piece(s)/minigrip bag |

Product characteristics

| | |
|----------------------------|----------------|
| Graduation | yes |
| Cap | flat cap |
| Closure type | push cap |
| Suitable for qPCR | yes |
| Base shape | conical |
| High or low profile | High Profile |
| Number of wells | 1 |
| Volume of work | 25 µl - 0,5 ml |
| Centrifugation max (RCF) | 8000 x g |
| Centrifugation temperature | 22 |
| Centrifugation time | 5 Min |
| Autoclavable (max. °C) | 121 °C |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
PCR single tube, 0.5 ml, PCR Performance
Tested, transparent, PP, flat cap

Page 2

Size

| | |
|-------------------|---------|
| Sample volume | 0.5 ml |
| Diameter | 7.9 mm |
| Width of product | 10 mm |
| Height of product | 29.5 mm |
| Length of product | 30.9 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Purity standard | PCR Performance Tested |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 2000 |
| Type of smallest subpackaging | minigrip bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 192000 |
| Depth of case | 348 mm |
| Width of case | 278 mm |
| Height of case | 140 mm |
| Case volume | 0.0135 cbm |
| Weight of product | 0.0004 kg |
| Weight of case | 1.1 kg |
| EAN of inner box | 4038917119986 |
| EAN case | 4038917016988 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT