




Product Specification

PCR plate without skirt, 96 well, white, Low Profile, 100 µl, PCR Performance Tested, PP

Page 1



Product description

Order number	72.1977.232
Product description	PCR plate without skirt, 96 well, white, Low Profile, 100 µl, PCR Performance Tested, material: PP, optimised for qPCR, 20 piece(s)/minigrip bag

Product characteristics

Closure type	film
Suitable for qPCR	yes
Base shape	conical
Construction of the plate edges	without skirt
High or low profile	Low Profile
Number of wells	96
Volume of work	5 µl - 0,1 ml
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
PCR plate without skirt, 96 well, white, Low Profile, 100 µl, PCR Performance Tested, PP

Page 2

Size

Sample volume	0.1 ml
Max. volume	100 µl
Width of product	77.7 mm
Length of product	117.7 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	white

Purity & certification

Product category	no medical device IVD
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	100
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	20
Piece (s) / outer case	100
Piece(s) / pallet	5000
Depth of case	137 mm
Width of case	137 mm
Height of case	493 mm
Case volume	0.0093 cbm
Weight of product	0.0106 kg
Weight of case	1.5 kg
EAN of inner box	4038917405133
EAN case	4038917405126

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT