




## Product Specification

Rack, PC, format: 12 x 4, suitable for tubes,  
square cuvettes, all S-Monovette® diameters

Page 1



### Product description

Order number	93.1431.001
Product description	Rack, material: PC, transparent, format: 12 x 4, (LxWxH): 257 x 90 x 40 mm, for 48 collection tubes, suitable for tubes, square cuvettes, all S-Monovette® diameters, 1 piece(s)/case

### Product characteristics

Suitable for	tubes, square cuvettes, all S-Monovette® diameters
Storage capacity	48
Format	12 x 4
Autoclavable (max. °C)	121 °C

### Size

Width of product	90 mm
Height of product	40 mm
Length of product	257 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-06-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Rack, PC, format: 12 x 4, suitable for tubes,  
square cuvettes, all S-Monovette® diameters**

**Material & colours**

Product material Polycarbonate (PC)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Batched no

**Packaging**

Minimum order qty. 1  
Type of smallest case  
subpackaging  
Piece(s) / inner box 1  
Piece (s) / outer case 40  
Piece(s) / pallet 1200  
Depth of box 92 mm  
Width of box 42 mm  
Height of box 259 mm  
Depth of case 398 mm  
Width of case 398 mm  
Height of case 340 mm  
Case volume 0.0539 cbm  
Weight of product 0.2036 kg  
Weight of case 9.88 kg  
EAN of inner box 4038917129107  
EAN case 4038917021968

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

