



### **Product description**

Order number	93.1707
Product description	Rack, material: PP, blue, format: 6 x 3, (LxWxH): 295 x 114 x 65 mm, for 18 collection tubes, suitable for 50 ml centrifuge tubes, 1 piece(s)/bag

### **Product characteristics**

Suitable for	50 ml centrifuge tubes
Storage capacity	18
Format	6 x 3
Autoclavable (max. °C)	121 °C

### **Size**

Width of product	114 mm
Height of product	65 mm
Length of product	295 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Rack, PP, format: 6 x 3, suitable for 50 ml**  
**centrifuge tubes**

**Page 2**

**Material & colours**

Product material Polypropylene (PP)  
Colour of product blue

**Purity & certification**

Product category no medical device | IVD  
Batched no

**Packaging**

Minimum order qty. 1  
Type of smallest subpackaging bag  
Piece(s) / inner box 1  
Piece (s) / outer case 10  
Piece(s) / pallet 540  
Depth of case 125 mm  
Width of case 75 mm  
Height of case 300 mm  
Case volume 0.0028 cbm  
Weight of product 0.2298 kg  
Weight of case 2.81 kg  
EAN of inner box 4038917347501  
EAN case 4038917347495

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**