



Product description

Order number	95.64.952
Product description	Storage box, slip-on lid, material: cardboard, white, format: 9 x 9, for 81 collection tubes, suitable for tubes Ø 13 mm

Product characteristics

Closure type	slip-on lid
Suitable for	tubes Ø 13 mm
Storage capacity	81
Format	9 x 9

Size

Width of product	135 mm
Height of product	105 mm
Length of product	135 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Storage box, slip-on lid, cardboard, format: 9 x 9,
for 81 collection tubes**

Page 2

Material & colours

Product material	cardboard
Colour of product	white
Closure material	cardboard
Colour of cap	white

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	60
Piece(s) / inner box	1
Piece (s) / outer case	60
Piece(s) / pallet	720
Case volume	0.1387 cbm
Weight of product	0.1298 kg
Weight of case	7.82 kg
EAN of inner box	4038917319638
EAN case	4038917571081

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT