



### **Product description**

Order number	95.64.964
Product description	Storage box, slip-on lid, material: cardboard, white, format: 8 x 8, for 64 collection tubes, suitable for tubes Ø 15 mm, 1 piece(s)/case

### **Product characteristics**

Label/ Print	with print
Closure type	slip-on lid
Suitable for	tubes Ø 15 mm
Storage capacity	64
Format	8 x 8

### **Size**

Width of product	135 mm
Height of product	45 mm
Length of product	135 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	cardboard
Colour of product	white
Closure material	cardboard
Colour of cap	white

### **Purity & certification**

Product category	no medical device   IVD
Batched	yes

### **Packaging**

Minimum order qty.	120
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	120
Piece(s) / pallet	1920
Case volume	0.1134 cbm
Weight of product	0.0908 kg
Weight of case	10.95 kg
EAN of inner box	4038917247412
EAN case	4038917247405

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

