



### **Product description**

Order number	65.729.301
Product description	Screw cap, red, suitable for tubes Ø 15.3 mm, 1,000 piece(s)/bag

### **Product characteristics**

Closure type	screw cap
Suitable for	tubes Ø 15.3 mm

### **Material & colours**

Product material	High Density Polyethylene (HD-PE)
Colour of product	red

### **Purity & certification**

Product category	no medical device   IVD
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	240000
Case volume	0.0309 cbm
Weight of product	0.001 kg
Weight of case	5.36 kg
EAN of inner box	4038917368339
EAN case	4038917368322

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

