



### **Product description**

Order number	65.791
Product description	Push cap, natural, suitable for tubes Ø 24.8 mm, 1,000 piece(s)/bag

### **Product characteristics**

Cap	standard cap
Closure type	push cap
Suitable for	tubes Ø 24.8 mm

### **Material & colours**

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	84000
Depth of case	388 mm
Width of case	378 mm
Height of case	235 mm
Case volume	0.0345 cbm
Weight of product	0.002 kg
Weight of case	4.61 kg
EAN of inner box	4038917084147
EAN case	4038917002998

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

