



### **Product description**

|                     |   |
|---------------------|---|
| Order number        | 65.802  |
| Product description | Push cap, natural, suitable for tubes Ø 15.7 mm, flat, 500 piece(s)/bag |

### **Product characteristics**

|              |                 |
|--------------|-----------------|
| Cap          | flat cap        |
| Closure type | push cap        |
| Suitable for | tubes Ø 15.7 mm |

### **Material & colours**

|                   |                                  |
|-------------------|----------------------------------|
| Product material  | Low Density Polyethylene (LD-PE) |
| Colour of product | natural                          |

### **Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | no                      |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 10000         |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 500           |
| Piece (s) / outer case        | 10000         |
| Piece(s) / pallet             | 360000        |
| Depth of case                 | 398 mm        |
| Width of case                 | 398 mm        |
| Height of case                | 300 mm        |
| Case volume                   | 0.0475 cbm    |
| Weight of product             | 0.0006 kg     |
| Weight of case                | 6.51 kg       |
| EAN of inner box              | 4038917084185 |
| EAN case                      | 4038917003049 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

