



Product description

Order number	65.802.100
Product description	Push cap, natural, suitable for tubes Ø 15.7 mm, flat, 500 piece(s)/bag

Product characteristics

Cap	flat cap
Closure type	push cap
Suitable for	tubes Ø 15.7 mm

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	natural

Purity & certification

Product category	no medical device IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	10000
Piece(s) / pallet	360000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0006 kg
Weight of case	6.71 kg
EAN of inner box	4038917989732
EAN case	4038917989749

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

