



### Product description

Order number	55.459
Product description	Tube, working volume: 13 ml, (LxØ): 100 x 16 mm, material: PS, round base, transparent, push cap, 500 piece(s)/bag

### Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	13 ml
Centrifugation max (RCF)	2500 x g

### Size

Diameter	16 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

### **Purity & certification**

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0046 kg
Weight of case	5.28 kg
EAN of inner box	4038917081672
EAN case	4038917000789

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

