



### **Product description**

Order number	55.461
Product description	Tube, working volume: 14 ml, (LxØ): 119 x 15.5 mm, material: PS, round base, transparent, push cap, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	14 ml
Centrifugation max (RCF)	3500 x g

### **Size**

Diameter	15.5 mm
Length of product	119 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material                      Polystyrene (PS)  
Colour of product                      transparent

### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      yes

### **Packaging**

Minimum order qty.                      1000  
Type of smallest  
subpackaging                              bag  
Piece(s) / inner box                      500  
Piece (s) / outer case                      1000  
Piece(s) / pallet                              30000  
Depth of case                              292 mm  
Width of case                                464 mm  
Height of case                               413 mm  
Case volume                                0.0559 cbm  
Weight of product                        0.0058 kg  
Weight of case                               6.29 kg  
EAN of inner box                            4038917081696  
EAN case                                      4038917000802

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

