



Product description

Order number	55.463
Product description	Tube, working volume: 14 ml, (LxØ): 105 x 16.8 mm, material: PS, round base, transparent, push cap, 250 piece(s)/StackPack

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	14 ml
Centrifugation max (RCF)	3500 x g

Size

Diameter	16.8 mm
Length of product	105 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Satisfies the requirement sealed at 0.6 bar/10 min
Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1000
Type of smallest
subpackaging StackPack
Piece(s) / inner box 250
Piece (s) / outer case 1000
Piece(s) / pallet 42000
Depth of box 360 mm
Width of box 107 mm
Height of box 184 mm
Depth of case 388 mm
Width of case 378 mm
Height of case 235 mm
Case volume 0.0345 cbm
Weight of product 0.0062 kg
Weight of case 6.95 kg
EAN of inner box 4038917081719
EAN case 4038917000826

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

