



Product description

Order number	55.466
Product description	Tube, working volume: 11.5 ml, (LxØ): 100 x 15.7 mm, material: PS, round base, transparent, push cap, 250 piece(s)/StackPack

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	11,5 ml
Centrifugation max (RCF)	3500 x g

Size

Diameter	15.7 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Satisfies the requirement sealed at 0.6 bar/10 min
Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1000
Type of smallest StackPack
subpackaging
Piece(s) / inner box 250
Piece (s) / outer case 1000
Piece(s) / pallet 48000
Depth of box 336 mm
Width of box 102 mm
Height of box 174 mm
Depth of case 388 mm
Width of case 378 mm
Height of case 215 mm
Case volume 0.0315 cbm
Weight of product 0.0055 kg
Weight of case 6.1 kg
EAN of inner box 4038917081740
EAN case 4038917000857

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

