



Product description

Order number	55.468
Product description	Tube, working volume: 13 ml, (LxØ): 95 x 16.8 mm, material: PS, round base, transparent, push cap, 250 piece(s)/StackPack

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	13 ml
Centrifugation max (RCF)	4000 x g

Size

Diameter	16.8 mm
Length of product	95 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Satisfies the requirement sealed at 0.6 bar/10 min
Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1000
Type of smallest
subpackaging StackPack
Piece(s) / inner box 250
Piece (s) / outer case 1000
Piece(s) / pallet 48000
Depth of box 366 mm
Width of box 97 mm
Height of box 186 mm
Depth of case 388 mm
Width of case 378 mm
Height of case 215 mm
Case volume 0.0315 cbm
Weight of product 0.0058 kg
Weight of case 6.38 kg
EAN of inner box 4038917081764
EAN case 4038917000871

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

