



Product description

Order number 55.475.030
Product description Tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PS, round base, transparent, 500 piece(s)/StackPack

Product characteristics

Graduation no
Base shape round base
Volume of work 5 ml
Centrifugation max (RCF) 4000 x g

Size

Diameter 13 mm
Length of product 75 mm

Material & colours

Product material Polystyrene (PS)
Colour of product transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	StackPack
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	112000
Depth of box	416 mm
Width of box	77 mm
Height of box	193 mm
Depth of case	428 mm
Width of case	383 mm
Height of case	175 mm
Case volume	0.0287 cbm
Weight of product	0.0031 kg
Weight of case	6.89 kg
EAN of inner box	4038917357883
EAN case	4038917357876

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

