



Product description

Order number	55.476.061
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PS, round base, transparent, push cap, suitable for Vitek and Vitek 2 from bioMérieux, 500 piece(s)/StackPack

Product characteristics

Graduation	no
Closure type	push cap
Suitable for	Vitek and Vitek 2 from bioMérieux
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	4000 x g

Size

Diameter	12 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Satisfies the requirement sealed at 0.6 bar/10 min
Product category In vitrodiagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 2000
Type of smallest StackPack
subpackaging
Piece(s) / inner box 500
Piece (s) / outer case 2000
Piece(s) / pallet 120000
Depth of box 387 mm
Width of box 77 mm
Height of box 176 mm
Depth of case 398 mm
Width of case 358 mm
Height of case 175 mm
Case volume 0.0249 cbm
Weight of product 0.0028 kg
Weight of case 6.2 kg
EAN of inner box 4038917107631
EAN case 4038917011839

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

