



Product description

Order number	55.516.177
Product description	Tube, working volume: 8 ml, (LxØ): 100 x 13 mm, material: PP, round base, brown, push cap, 500 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	8 ml
Centrifugation max (RCF)	3500 x g
Autoclavable (max. °C)	121 °C

Size

Diameter	13 mm
Height of product	100 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product brown

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 500
Type of smallest subpackaging bag
Piece(s) / inner box 500
Piece (s) / outer case 2000
Piece(s) / pallet 56000
Case volume 0.0649 cbm
Weight of product 0.0034 kg
Weight of case 7.85 kg
EAN of inner box 4038917108980
EAN case 4038917031981

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

