



### **Product description**

Order number	55.525.505
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 13 mm, material: PP, round base, transparent, push cap, 1,000 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g

### **Size**

Diameter	13 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

### **Purity & certification**

Satisfies the requirement sealed at 0.6 bar/10 min  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

### **Packaging**

Minimum order qty. 2000  
Type of smallest subpackaging bag  
Piece(s) / inner box 1000  
Piece (s) / outer case 2000  
Piece(s) / pallet 84000  
Depth of case 398 mm  
Width of case 398 mm  
Height of case 250 mm  
Case volume 0.0396 cbm  
Weight of product 0.0026 kg  
Weight of case 5.78 kg  
EAN of inner box 4038917375955  
EAN case 4038917375948

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

