



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 55.526.046   |
| Product description | Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, brown, push cap, with graduations, 1,000 piece(s)/bag |

### **Product characteristics**

|                          |            |
|--------------------------|------------|
| Graduation               | yes        |
| Closure type             | push cap   |
| Base shape               | round base |
| Volume of work           | 5 ml       |
| Centrifugation max (RCF) | 5000 x g   |
| Autoclavable (max. °C)   | 121 °C     |

### **Size**

|                   |       |
|-------------------|-------|
| Diameter          | 12 mm |
| Height of product | 75 mm |
| Length of product | 75 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product brown

### **Purity & certification**

Product category no medical device | IVD  
Batched yes

### **Packaging**

Minimum order qty. 4000  
Type of smallest subpackaging bag  
Piece(s) / inner box 1000  
Piece (s) / outer case 4000  
Piece(s) / pallet 64000  
Case volume 0.0813 cbm  
Weight of product 0.0021 kg  
Weight of case 10.34 kg  
EAN of inner box 4038917109642  
EAN case 4038917068581

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

