



Product description

| | |
|---------------------|---|
| Order number | 55.526.521 |
| Product description | Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, blue, push cap, 500 piece(s)/StackPack |

Product characteristics

| | |
|--------------------------|------------|
| Graduation | no |
| Closure type | push cap |
| Base shape | round base |
| Volume of work | 5 ml |
| Centrifugation max (RCF) | 5000 x g |
| Autoclavable (max. °C) | 121 °C |

Size

| | |
|-------------------|-------|
| Diameter | 12 mm |
| Height of product | 75 mm |
| Length of product | 75 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product blue

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 500
Type of smallest subpackaging StackPack
Piece(s) / inner box 500
Piece (s) / outer case 4000
Piece(s) / pallet 120000
Case volume 0.046 cbm
Weight of product 0.0021 kg
Weight of case 10.25 kg
EAN of inner box 4038917109758
EAN case 4038917031950

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

