



### **Product description**

Order number 55.526.532  
Product description Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, yellow, push cap, 1,000 piece(s)/bag

### **Product characteristics**

Graduation no  
Closure type push cap  
Base shape round base  
Volume of work 5 ml  
Centrifugation max (RCF) 5000 x g  
Autoclavable (max. °C) 121 °C

### **Size**

Diameter 12 mm  
Height of product 75 mm  
Length of product 75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product yellow

### **Purity & certification**

Product category no medical device | IVD  
Batched yes

### **Packaging**

Minimum order qty. 1000  
Type of smallest subpackaging bag  
Piece(s) / inner box 1000  
Piece (s) / outer case 4000  
Piece(s) / pallet 64000  
Case volume 0.0813 cbm  
Weight of product 0.0021 kg  
Weight of case 10.3 kg  
EAN of inner box 4038917109826  
EAN case 4038917032247

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

