



Product description

Order number	57.462.300
Product description	Tube, working volume: 12 ml, (LxØ): 110 x 17 mm, material: PS, conical base, transparent, push cap, with graduations, 500 piece(s)/bag

Product characteristics

Graduation	yes
Closure type	push cap
Base shape	conical base
Volume of work	12 ml
Centrifugation max (RCF)	2500 x g

Size

Diameter	17 mm
Height of product	110 mm
Length of product	110 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	32000
Case volume	0.0558 cbm
Weight of product	0.0061 kg
Weight of case	6.61 kg
EAN of inner box	4038917367400
EAN case	4038917367394

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

