



Product description

Order number 58.485
Product description Tube, working volume: 7 ml, (LxØ): 50 x 16 mm, material: PS, flat base, transparent, push cap, 1,000 piece(s)/bag

Product characteristics

Graduation no
Closure type push cap
Base shape flat base
Volume of work 7 ml
Centrifugation max (RCF) 2500 x g

Size

Diameter 16 mm
Length of product 50 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|------------------|
| Product material | Polystyrene (PS) |
| Colour of product | transparent |

Purity & certification

| | |
|---------------------------|--------------------------|
| Satisfies the requirement | sealed at 0.6 bar/10 min |
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 2000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1000 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 72000 |
| Depth of case | 398 mm |
| Width of case | 398 mm |
| Height of case | 300 mm |
| Case volume | 0.0475 cbm |
| Weight of product | 0.0028 kg |
| Weight of case | 6.34 kg |
| EAN of inner box | 4038917081924 |
| EAN case | 4038917001038 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

