



Product description

Order number	58.491.300
Product description	Tube, working volume: 12 ml, (LxØ): 95 x 16.5 mm, material: PS, flat base, transparent, push cap, 500 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	flat base
Volume of work	12 ml

Size

Diameter	16.5 mm
Height of product	95 mm
Length of product	95 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest
subpackaging bag
Piece(s) / inner box 500
Piece (s) / outer case 1000
Piece(s) / pallet 36000
Case volume 0.0492 cbm
Weight of product 0.0056 kg
Weight of case 6.07 kg
EAN of inner box 4038917367493
EAN case 4038917367486

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

