



Product description

Order number	58.505
Product description	Tube, working volume: 12 ml, (LxØ): 95 x 16.5 mm, material: PP, flat base, transparent, push cap, 500 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	flat base
Volume of work	12 ml
Centrifugation max (RCF)	3500 x g

Size

Diameter	16.5 mm
Length of product	95 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Satisfies the requirement sealed at 0.6 bar/10 min
Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1000
Type of smallest bag
subpackaging
Piece(s) / inner box 500
Piece (s) / outer case 1000
Piece(s) / pallet 36000
Depth of case 398 mm
Width of case 398 mm
Height of case 265 mm
Case volume 0.042 cbm
Weight of product 0.0048 kg
Weight of case 5.38 kg
EAN of inner box 4038917082020
EAN case 4038917001137

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

