



Product description

| | |
|---------------------|---|
| Order number | 58.505.046 |
| Product description | Tube, working volume: 12 ml, (LxØ): 95 x 16.5 mm, material: PP, flat base, brown, push cap, brown, cap enclosed, 500 piece(s)/bag |

Product characteristics

| | |
|----------------|-----------|
| Graduation | no |
| Closure type | push cap |
| Base shape | flat base |
| Volume of work | 12 ml |

Size

| | |
|-------------------|---------|
| Diameter | 16.5 mm |
| Length of product | 95 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | brown |
| Closure material | Low Density Polyethylene (LD-PE) |
| Colour of cap | brown |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 32000 |
| Case volume | 0.0558 cbm |
| Weight of product | 0.0064 kg |
| Weight of case | 6.89 kg |
| EAN of inner box | 4038917225557 |
| EAN case | 4038917225540 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

