



### **Product description**

Order number 58.536.300  
Product description Tube, working volume: 7 ml, (LxØ): 50 x 16 mm, material: PP, flat base, transparent, push cap, 1,000 piece(s)/bag

### **Product characteristics**

Graduation no  
Closure type push cap  
Base shape flat base  
Volume of work 7 ml  
Autoclavable (max. °C) 121 °C

### **Size**

Diameter 16 mm  
Height of product 50 mm  
Length of product 50 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

### **Purity & certification**

Product category	no medical device   IVD
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	72000
Case volume	0.0492 cbm
Weight of product	0.0024 kg
Weight of case	5.24 kg
EAN of inner box	4038917367554
EAN case	4038917367547

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

