



### **Product description**

Order number	58.618.008
Product description	Tube, working volume: 12 ml, (LxØ): 95 x 16.5 mm, material: PS, skirted conical base, brown, push cap, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	skirted conical base
Volume of work	12 ml

### **Size**

Diameter	16.5 mm
Height of product	95 mm
Length of product	95 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material                      Polystyrene (PS)  
Colour of product                      brown

### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      yes

### **Packaging**

Minimum order qty.                      1000  
Type of smallest  
subpackaging                              bag  
Piece(s) / inner box                      500  
Piece (s) / outer case                      1000  
Piece(s) / pallet                              36000  
Case volume                                0.0492 cbm  
Weight of product                        0.0056 kg  
Weight of case                              6.1 kg  
EAN of inner box                            4038917115148  
EAN case                                      4038917032414

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

