



Product description

Order number	58.619
Product description	Tube, working volume: 12 ml, (LxØ): 95 x 16.5 mm, material: PP, skirted conical base, transparent, push cap, 500 piece(s)/bag

Product characteristics

Graduation	no
Closure type	push cap
Base shape	skirted conical base
Volume of work	12 ml
Autoclavable (max. °C)	121 °C

Size

Diameter	16.5 mm
Length of product	95 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest subpackaging bag
Piece(s) / inner box 500
Piece (s) / outer case 1000
Piece(s) / pallet 36000
Case volume 0.0492 cbm
Weight of product 0.0056 kg
Weight of case 6.1 kg
EAN of inner box 4038917082945
EAN case 4038917032421

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

