



**Product Specification**  
**Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, without cap, 1,000 piece(s)/bag**

**Page 1**



**Product description**

Order number	60.504.031
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 75 x 13 mm, round base, transparent, material: PP, without cap, 1,000 piece(s)/bag, 2,000 piece(s)/carton

**Product characteristics**

Graduation	no
Closure type	screw cap
Base shape	round base
Volume of work	5 ml

**Size**

Diameter	13 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 5 ml, (LxØ): 75 x 13 mm, round base, PP, without cap, 1,000 piece(s)/bag**

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Batched yes

**Packaging**

Minimum order qty. 1000  
Type of smallest subpackaging bag  
Piece(s) / inner box 1000  
Piece (s) / outer case 2000  
Piece(s) / pallet 72000  
Depth of case 458 mm  
Width of case 408 mm  
Height of case 260 mm  
Case volume 0.0486 cbm  
Weight of product 0.0026 kg  
Weight of case 5.75 kg  
EAN of inner box 4038917412735  
EAN case 4038917412728

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

