



Product description

| | |
|---------------------|---|
| Order number | 60.506.024 |
| Product description | Tube, working volume: 10 ml, (LxØ): 79 x 16 mm, material: PP, round base, transparent, screw cap, without cap, 500 piece(s)/bag |

Product characteristics

| | |
|--------------------------|------------|
| Graduation | no |
| Closure type | screw cap |
| Base shape | round base |
| Volume of work | 10 ml |
| Centrifugation max (RCF) | 5000 x g |
| Autoclavable (max. °C) | 121 °C |

Size

| | |
|-------------------|-------|
| Diameter | 16 mm |
| Length of product | 79 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 500
Type of smallest subpackaging bag
Piece(s) / inner box 500
Piece (s) / outer case 2000
Piece(s) / pallet 32000
Case volume 0.0813 cbm
Weight of product 0.0032 kg
Weight of case 7.39 kg
EAN of inner box 4038917108713
EAN case 4038917035873

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

