



### Product description

Order number	60.540.686
Product description	Screw cap tube, working volume: 13 ml, (LxØ): 101 x 16.5 mm, material: PP, round base with skirted base, transparent, screw cap, natural, cap assembled, sterile, 500 piece(s)/bag

### Product characteristics

Graduation	no
Closure type	screw cap
Base shape	round base with skirted base
Volume of work	13 ml
Centrifugation max (RCF)	4000 x g

### Size

Diameter	16.5 mm
Length of product	101 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

### **Purity & certification**

Satisfies the requirement	IATA
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### **Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	500
Piece(s) / pallet	22000
Depth of case	568 mm
Width of case	348 mm
Height of case	155 mm
Case volume	0.0306 cbm
Weight of product	0.0066 kg
Weight of case	3.8 kg
EAN of inner box	4038917110501
EAN case	4038917012898

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

