



Product description

Order number	60.542.307
Product description	Screw cap tube, working volume: 8 ml, (LxØ): 57 x 16.5 mm, material: PP, flat base, transparent, screw cap, natural, cap enclosed, 500 piece(s)/bag

Product characteristics

Graduation	no
Closure type	screw cap
Base shape	flat base
Volume of work	8 ml
Centrifugation max (RCF)	4000 x g
Autoclavable (max. °C)	121 °C

Size

Diameter	16.5 mm
Length of product	57 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	natural
O-ring material	Ethylene propylene diene rubber (EPDM)

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	40000
Case volume	0.0382 cbm
Weight of product	0.0043 kg
Weight of case	4.81 kg
EAN of inner box	4038917376044
EAN case	4038917376037

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

