



### **Product description**

Order number	60.545
Product description	Screw cap tube, working volume: 30 ml, (LxØ): 107 x 25 mm, material: PP, conical base, transparent, screw cap, natural, cap assembled, sterile, 50 piece(s)/bag

### **Product characteristics**

Graduation	no
Closure type	screw cap
Base shape	conical base
Volume of work	30 ml
Centrifugation max (RCF)	8800 x g

### **Size**

Diameter	25 mm
Length of product	107 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

### **Purity & certification**

Satisfies the requirement	ADR
Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### **Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	14000
Depth of case	588 mm
Width of case	398 mm
Height of case	255 mm
Case volume	0.0597 cbm
Weight of product	0.0113 kg
Weight of case	6.31 kg
EAN of inner box	4038917082389
EAN case	4038917001489

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

