



Product description

Order number	60.551.023
Product description	Screw cap tube, working volume: 10 ml, (LxØ): 79 x 16 mm, material: PP, round base with skirted base, transparent, screw cap, without cap, 500 piece(s)/bag

Product characteristics

Graduation	no
Closure type	screw cap
Base shape	round base with skirted base
Volume of work	10 ml
Centrifugation max (RCF)	4000 x g
Autoclavable (max. °C)	121 °C

Size

Diameter	16 mm
Height of product	79 mm
Length of product	79 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 2000
Type of smallest subpackaging bag
Piece(s) / inner box 500
Piece (s) / outer case 2000
Piece(s) / pallet 32000
Case volume 0.0813 cbm
Weight of product 0.0024 kg
Weight of case 5.9 kg
EAN of inner box 4038917317801
EAN case 4038917317795

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

