



Product description

Order number	60.558
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 57 x 15.3 mm, material: PP, skirted conical base, transparent, screw cap, natural, cap enclosed, 1,000 piece(s)/bag

Product characteristics

Graduation	no
Closure type	screw cap
Base shape	skirted conical base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g

Size

Diameter	15.3 mm
Length of product	57 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

Purity & certification

Satisfies the requirement	IATA
Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	1000
Piece(s) / pallet	48000
Depth of case	388 mm
Width of case	378 mm
Height of case	215 mm
Case volume	0.0315 cbm
Weight of product	0.0037 kg
Weight of case	4.14 kg
EAN of inner box	4038917082464
EAN case	4038917001540

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

