



Product description

Order number	60.596.001
Product description	Screw cap tube, working volume: 60 ml, (LxØ): 126 x 30 mm, material: PP, flat base, transparent, screw cap, natural, cap assembled, sterile, 25 piece(s)/bag

Product characteristics

Graduation	no
Cap	standard screw cap
Closure type	screw cap
Base shape	flat base
Volume of work	60 ml
Autoclavable (max. °C)	121 °C

Size

Diameter	30 mm
Length of product	126 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	natural

Purity & certification

Satisfies the requirement	sealed at 0.6 bar/10 min
Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	250
Type of smallest subpackaging	bag
Piece(s) / inner box	25
Piece (s) / outer case	250
Piece(s) / pallet	7000
Depth of case	588 mm
Width of case	398 mm
Height of case	255 mm
Case volume	0.0597 cbm
Weight of product	0.015 kg
Weight of case	4.44 kg
EAN of inner box	4038917114134
EAN case	4038917014199

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

