



Product description

Order number	60.610.100
Product description	Screw cap tube, working volume: 10 ml, (LxØ): 92 x 15.3 mm, material: PP, round base, transparent, screw cap, transparent, cap assembled, 100 piece(s)/bag

Product characteristics

Graduation	no
Closure type	screw cap
Base shape	round base
Volume of work	10 ml
Centrifugation max (RCF)	3000 x g

Size

Diameter	15.3 mm
Length of product	92 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	transparent
Tube material	Polypropylene (PP)

Purity & certification

Satisfies the requirement	IATA
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	36000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0052 kg
Weight of case	5.92 kg
EAN of inner box	4038917114776
EAN case	4038917014380

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

