



### **Product description**

Order number	60.9921.524
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 50 x 16 mm, material: PP, flat base, translucent, screw cap, yellow, cap assembled, sterile, 250 piece(s)/bag

### **Product characteristics**

Graduation	no
Cap	standard screw cap
Closure type	screw cap
Base shape	flat base
Volume of work	5 ml
Centrifugation max (RCF)	5500 x g
Centrifugation temperature	21
Centrifugation time	15 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### Size

Diameter	16 mm
Length of product	50 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	translucent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Tube material	Polypropylene (PP)

### Purity & certification

Satisfies the requirement	sealed at 18"hg/10min
Latex-free	yes
Purity standard	sterile
Pyrogen-free	no
Safety data sheet	no

### Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	2000
Piece(s) / pallet	96000
Depth of case	608 mm
Width of case	264 mm
Height of case	324 mm
Case volume	0.0519 cbm
Weight of case	7.4 kg
EAN of inner box	4038917150804
EAN case	4038917037150

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

