



### **Product description**

Order number	60.9922.934
Product description	Screw cap tube, working volume: 12 ml, (LxØ): 99 x 16 mm, material: PS, round base, transparent, screw cap, yellow, cap assembled, sterile, 100 piece(s)/bag

### **Product characteristics**

Graduation	no
Cap	standard screw cap
Closure type	screw cap
Base shape	round base
Volume of work	12 ml
Centrifugation max (RCF)	5000 x g
Centrifugation temperature	21
Centrifugation time	15 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### Size

Diameter	16 mm
Length of product	99 mm

### Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Tube material	Polystyrene (PS)

### Purity & certification

Satisfies the requirement	sealed at 18"hg/10min
Latex-free	yes
Purity standard	sterile
Pyrogen-free	no
Safety data sheet	no

### Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	40000
Depth of case	604 mm
Width of case	264 mm
Height of case	373 mm
Case volume	0.06 cbm
Weight of case	7.3 kg
EAN of inner box	4038917153416
EAN case	4038917037402

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

