



### **Product description**

Order number	62.476.019
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PS, round base, transparent, push cap, with print, label/print: orange, with graduations, 1,000 piece(s)/bag

### **Product characteristics**

Label/ Print	with print
Colour of print/label	orange
Graduation	yes
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	3500 x g

### **Size**

Diameter	12 mm
Height of product	75 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material                      Polystyrene (PS)  
Colour of product                      transparent

### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      yes

### **Packaging**

Minimum order qty.                      1000  
Type of smallest  
subpackaging                              bag  
Piece(s) / inner box                      1000  
Piece (s) / outer case                      2000  
Piece(s) / pallet                              80000  
Case volume                                0.0382 cbm  
Weight of product                        0.0026 kg  
Weight of case                              6.2 kg  
EAN of inner box                            4038917107372  
EAN case                                      4038917032643

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

