



Product description

Order number	62.476.022
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PS, round base, transparent, push cap, with print, label/print: white, with graduations, 1,000 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	4000 x g

Size

Diameter	12 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polystyrene (PS)
Colour of product transparent

Purity & certification

Satisfies the requirement sealed at 0.6 bar/10 min
Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 2000
Type of smallest bag
subpackaging
Piece(s) / inner box 1000
Piece (s) / outer case 2000
Piece(s) / pallet 84000
Depth of case 398 mm
Width of case 398 mm
Height of case 250 mm
Case volume 0.0396 cbm
Weight of product 0.0028 kg
Weight of case 6.22 kg
EAN of inner box 4038917107402
EAN case 4038917011792

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

