



Product Specification
**Screw cap tube, 8 ml, (LxØ): 57 x 16.5 mm, PP,
with print**

Page 1



Product description

Order number	62.542.009
Product description	Screw cap tube, working volume: 8 ml, (LxØ): 57 x 16.5 mm, material: PP, flat base, transparent, screw cap, natural, cap assembled, with print, label/print: white, with graduations, 1,000 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	screw cap
Base shape	flat base
Volume of work	8 ml
Centrifugation max (RCF)	4000 x g
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**Screw cap tube, 8 ml, (LxØ): 57 x 16.5 mm, PP,
with print**

Page 2

Size

Diameter 16.5 mm
Length of product 57 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap natural

Purity & certification

Satisfies the requirement IATA
Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1000
Type of smallest subpackaging bag
Piece(s) / inner box 1000
Piece (s) / outer case 1000
Piece(s) / pallet 48000
Case volume 0.0309 cbm
Weight of product 0.0045 kg
Weight of case 4.99 kg
EAN of inner box 4038917111126
EAN case 4038917072236

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT