



### **Product description**

Order number	62.542.026
Product description	Screw cap tube, working volume: 8 ml, (LxØ): 57 x 16.5 mm, material: PP, flat base, transparent, screw cap, without cap, with print, label/print: white, with graduations, 1,000 piece(s)/bag

### **Product characteristics**

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	screw cap
Base shape	flat base
Volume of work	8 ml
Centrifugation max (RCF)	4000 x g
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 8 ml, (LxØ): 57 x 16.5 mm, PP,  
with print**

**Page 2**

**Size**

Diameter	16.5 mm
Height of product	57 mm
Length of product	57 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	1000
Piece(s) / pallet	60000
Case volume	0.0282 cbm
Weight of product	0.0046 kg
Weight of case	5.15 kg
EAN of inner box	4038917111218
EAN case	4038917033206

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**